

Indicazioni d'uso Aditus Flexible graft

L'innesto osseo flessibile di origine biominerale.

	Post extraction socket grafting	Post extraction socket grafting (No implant scheduled)	Peri implant (Simultaneous implant & graft)
Aditus flexible graft	✗	✗	★ ★
GBR for delayed implant placement	Sinus lift (Crestal Approach)	Sinus lift lateral approach	Ridge agumentation (Only graft for width)
★	✗	✗	★ ★ ★
Ridge agumentation (Khoury technique)	Ridge agumentation (Ring technique)	Bone cyst grafting	Block graft (Autogenous)
✗	✗	★	★

Legenda:

✗ = Non adatto ★ = adatto ★ ★ = molto adatto ★ ★ ★ = Perfetto